

The right EPA-registered surface disinfectant makes all the difference in the fight against coronavirus.

Cetylcide II[®]

Hospital-Grade
Disinfectant Concentrate

JUST ADD WATER



Powerful. Effective. Simple.

Follow CDC guidelines and use only hospital-grade surface disinfectants that are EPA-registered for use against viruses similar to SARS-CoV-2, the virus that causes COVID-19.



Kills 140+ pathogens

on surfaces including human coronavirus*



EPA Registered

meets EPA's Emerging Viral Pathogen Policy



One step to clean & disinfect



No rinsing required

Protect the people you care about. Kill germs on the surfaces they touch with Cetylcide-II® Disinfectant Concentrate.

Effective against 140+ pathogens on surfaces including:

- HIV
- MRSA
- H1N1
- E. coli
- Zika Virus
- Rotavirus
- Human Coronavirus*
- Porcine Parvovirus
- Norovirus
- Streptococcus Pneumoniae
- Bordetella
- Salmonella
- Influenza
- Hepatitis B
- Hepatitis C
- Measles
- Poliovirus Type 1

***COVID-19 is caused by SARS-CoV-2.**

Cetylcide II® has demonstrated effectiveness against viruses similar to SARS-CoV-2 on hard, non-porous surfaces. Therefore, Cetylcide II® can be used against SARS-CoV-2 when used in accordance with the directions for use against Rotavirus on hard, non-porous surfaces. Refer to the CDC website at <http://www.cdc.gov/coronavirus/2019-nCoV/index.html> for additional information.

SPRAY SURFACE AND LET DRY



"This disinfectant works great and smells great. I bought it because I can't find wipes anywhere and now I found something that works even better. I can clean my kitchen, bath, counters and floors all with this product and know that my family is safe. I'm buying more ASAP."

Simple to use, just add water and Cetylcide-II® is ready to use immediately











Cetylcide-II Disinfectant Concentrate is easily diluted with tap water. It's pre-measured for accurate mixing, with no waiting and no activation time. Once combined with tap water, the solution is immediately active and ready to use.

Use 1:64 dilution for most applications

Ounces of Product	Amount of Water
1/2 ounce	1 quart
1 ounce	1/2 gallon
2 ounces	1 gallon

Cetylcide-II Disinfectant Concentrate, a third-generation quaternary ammonium formula, provides broad-spectrum effectiveness without harsh alcohol, damaging bleach or corrosive hydrogen peroxide. It's used routinely and effectively in professional settings like hospitals, medical facilities, labs, schools, and public spaces. It's also safe for use at home.

Keep a wide variety of surfaces safe and clean, including:

-  Glass, Metal, Stainless Steel, Granite, Marble, Vinyl
-  Enameled Surfaces, Plastic Upholstery
-  Plastic/Plexiglass
-  Tables, Chairs, Floors
-  Bed and Bath Linens
-  Non-Food Contact Surfaces
-  Kitchen Countertops
-  Bathtubs/Showers
-  Desks
-  Athletic Equipment



Mixing Instructions

- 1 Water**
Fill Quart Spray Bottle with 1 quart of tap water.



- 2 + Concentrate**
Add 1/2 oz. (15mL) of Cetylcide-II Concentrate.



- 3 + Spray Nozzle**
Secure trigger spray to bottle.



- 4 Ready to Use**
Immediately, no activation period required.



Directions For Use

After diluting, use Cetylcide-II for cleaning and disinfecting.


To clean: Apply diluted solution to hard, non-porous surfaces with a sprayer or by thoroughly wetting surfaces with a cloth, mop, or sponge. For sprayer applications, use a coarse spray device. Spray 6-8 inches from the surface and rub with brush, sponge, or cloth. Store solution in closed container for up to 64 days.

To disinfect: Wet surfaces with diluted solution. Treated surfaces must remain wet for 10 minutes. Let air dry. No rinsing needed for non-food surfaces. Refer to product insert for complete instructions and safety precautions.

For complete Instructions-for-Use, visit the "Resources" page at www.cetyllite.com/dental

**Use Cetylcode-II®
Disinfectant Concentrate in:**

- Healthcare/Dental/
Veterinary Settings
- EMS Facilities/Vehicles
- Daycares/Nursing Homes
- Public Restrooms/
Rest Areas
- Health Clubs/Spas
- Residential
- Laundry
- Kitchens
- Bathrooms
- Fogging/Wet Misting



Germ's can't hide. Here's how Cetylcode-II® works.

This is the simple and effective solution for a broad spectrum of challenges. Positively charged Cetylcode-II is inherently drawn to surfaces, disinfecting and killing harmful microorganisms everywhere it reaches.

- Disinfectant
- Cleaner
- Detergent
- Deodorizer
- Virucide
- Sanitizer
- Fungicide
- Mildewstat

WHAT CUSTOMERS ARE SAYING:



"I feel safe from nasty germs with Cetylcode-II on the job!"



"I have been using this in my business. It's worth the cost."



"Excellent for cleaning. Works well. Would order again."



"Use it for work and it's what works the best for us."

Available In Multiple Sizes:



(0152)

32 oz concentrate
+ quart sprayer
Yields 64 quarts
(16 gallons)



(0153)

1/2 oz concentrate
+ quart sprayer
Yields 1 quart



(0158)

1/2 oz concentrate
(6/pk refills)
Yields 6 quarts



(0159)

1/2 oz concentrate
(12/pk refills)
Yields 12 quarts



Available for purchase
through medical and dental
suppliers and Amazon.com.



Cetylcode-II EPA Registration No. 61178-1-3150 | Cetylcode-II EPA Registration No. 3150-NJ-1
9051 River Road, Pennsauken, NJ 08110 | 856-665-6111 | www.cetylcode.com

Rev: CII 092020 | Printed in the USA



Proudly Made
in the USA