Abdruckspritze / Impression syringe / Seringue d'empreinte / Jeringa de impresión / Siringa per impronta / Seringa de impressão / Шприц для снятия слепка / Afdrukspuit / Strzykawka do mas wyciskowych / 注射器	REF	532, 533, 534
Dosier Messöffel / Dosing spoon / Cuillère-mesure de dosage / Cuchara de dosificación / Cucchiaino dosatore / Colher de dosagem / Дозировочная ложка / Doseerlepel / Łyżka dozująca / 量匙	REF	4610
Beflüttete Abdruckpads / Vented impression pads / Tampons d'empreinte aérés / Plantillas de impresión con ventilación / Tamponi per impronta ventilati / Tampões para impressão ventilados / Вентилируемые накладки для снятия слепка / Geventileerde afdrukpads / Wentylowane stożki ochronne do odcisków / 透气式取模垫	REF	5671, 5672, 5673, 5674, 5675
Abdruckpads / Impression pads / Tampons d'empreinte / Plantillas de impresión / Tamponi per impronta / Tampões para impressão / Накладки для снятия слепка / Afdrukpads / Stożki ochronne do odcisków / 取模垫	REF	565, 566

Technische Daten / Technical data / Données techniques / Datos técnicos / Dati tecnici / Dados técnicos / Технические данные / Technische gegevens / Dane techniczne / 技术资料

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(EN)

Caution!
The impression of the outer ear canal may only be taken with impression material by persons with appropriate qualifications. The processing information and precautions which are set out in these directions for use must be observed under all circumstances. Failure to observe these may result in irreparable damage to the ear or the eardrum. If damage or undesirable side effects occur despite the warnings, we request that you report this with a detailed description of your procedure and the damage to or side effects for the patient.

Intended use
Addition vulcanizing ear impression silicone.

Indication for use
Ear impression material for reproduction of the outer ear canal of a patient with hearing loss, Tinnitus aurium or impaired hearing for the provision and adjustment of a hearing device.

Contraindication
Side effects
Otoform® may not be used if the patient is allergic to one of the ingredients. There are no side-effects to be expected from this medical device if processed and used properly. However, immun reactions (e.g. allergies) or local disorders of sensation (e.g. irritation in the auditory canal) cannot be ruled out in principle. If you become aware of adverse side-effects – and in case of doubt – please be sure to advise us in any case and describe the circumstances and symptoms as accurately as possible. We will investigate any indication.

Patient groups
Otoform® must not be used for patients suffering from a perforated eardrum or injuries to the ear canal.

Target Patient Group
Patients who need ear care.

Intended users
ENT specialist, hearing aid audiologist.

Product description
The ear impression material consists of addition curing silicone with a soft, thixotropic consistency for impression technology which creates slight pressure. The material is practically shrink-proof and is physiologically well tolerated thanks to the platinum hardener system. Slight generation of heat can be felt while Otoform® is curing.

Processing
The ear canal and the eardrum must be examined before the impression is taken. No impression may be taken in the event of abnormal findings such as inflammation or perforation of the eardrum. Cleaning of the ear canal and removal of hairs by an ENT doctor is required in the event of significant earwax deposits and hairs in the ear canal. An impression pad must be placed at the end of the ear canal before the eardrum before the impression is taken in order to protect the eardrum. The two components of the impression material are removed from the containers in equal parts with measuring spoons and pressed onto a mixing tray. The two components are now rapidly kneaded together with a spatula until the mixture is a uniform colour. For hygiene reasons, kneading by hand is only recommended with gloves (no latex gloves!). The kneaded mass is now filled into an impression syringe. The impression material is slowly applied into the ear canal and the concha. Please take care to ensure that the tip of the impression syringe is always in direct contact with the impression mass during application. Once the material has hardened, the impression can be carefully removed from the ear. The ear canal must be examined again after the impression has been taken. Carefully close the tins after use.

Disinfection
The Otoform® impressions can be disinfected with OtoVita® Professional disinfectant spray in accordance with the information in the associated fact sheet.

- Caution**
- The product is exclusively intended for the above-mentioned range of use and shall only be used by professionally qualified and instructed persons.
 - Wear protective clothing! Set material is chemically resistant. Marks on your clothing cannot be removed.
 - Carefully close the tins after use.
 - The chemical reaction between the two components may be inhibited by earwax and the use of latex gloves, hand lotions, cleaning agents, etc. For this reason, we recommend that you wear commercially available polyethylene gloves.
 - The setting times can be increased by lower temperatures or reduced by higher temperatures.

Serious incidents
All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

Single-use advice
After taking an impression, the hardened ear impression silicone cannot be used for taking another impression owing to the physical properties.

Disposal
Disposal of contents / container in accordance with local / regional / national / international regulations and safety data sheet.

Hazard identification
Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

Storage
Storage at room temperature 18°C–28°C (64°F–82°F).

Lot number / Expiry date
The lot number and the expiry date are indicated on the external packaging as well as on each tub and bucket. Both components A and B of a packaging must only be used in this combination. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

See the last page for technical data and required specifications.
Dreve Otoplastik GmbH is not liable for any damages caused by improper application of the product.

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