



RECHARGEABLE

PATENTED

Lighted Ear Cures and Tweezers with Magnification



30 MINUTES

Charging time

50 MINUTES

Battery life when fully charged



The LightBeam® ear curettes system is intended to aid in the removal of cerumen and foreign objects from the ear canal.

The LightBeam complete kit includes a rechargeable base for the light source, eight different types of curette tips (six each), two tweezers, and a magnifier that fits on the curette. These instruments are all made from clear polycarbonate plastic. The instruments are disposable and intended for single use. All parts are numbered for easy reordering.

• Easily rechargeable light source

reduce light source replacement costs

• Patented threaded curette connection

no concerns over accidental ejection

• Innovative and patented tweezer attachment

• Eight different curette tips

• Competitively priced

LightBeam Ear Cures

1-800-328-4757 | info@warnertechcare.com



Another quality product by
EAR TECHNOLOGY CORPORATION®

Ear Technology Corporation's Johnson City facilities have a Quality Management System certified to ISO 9001: 2015 for the Design, Manufacture, Sale, and Distribution of Hearing Healthcare Products.
Certified—Service Disabled Veteran Owned Small Business

PATENTED

Lighted Ear Curettes and Tweezers with Magnification

lightbeamcurettes.com

RECHARGEABLE



#1 | LIGHTED OVAL LOOP



#2 | LIGHTED SMALL SPOON



#3 | LIGHTED SCOOP



#4 | LIGHTED SMALL TRI-LOOP



#5 | LIGHTED LARGE CIRCLE LOOP



#6 | LIGHTED LARGE SPOON



#7 | LIGHTED SMALL CIRCLE LOOP



#8 | LIGHTED NARROW SPOON



#9 | LIGHTED TWEEZERS

COMPLETE STARTER KIT

626000

50* total curettes + 1 Light Source + 1 Magnifier
+ Charging Base, USB Cable, Adapter

* Six (6) each of 8 different curettes + 2 Tweezers

REFILL KITS

Part Number

50 pk	#1 Lighted Oval Loop	626001
50 pk	#2 Lighted Small Spoon	626002
50pk	#3 Lighted Scoop	626003
50 pk	#4 Lighted Small Tri-Loop	626004
50 pk	#5 Lighted Large Circle Loop	626005
50pk	#6 Lighted Large Spoon	626006
50 pk	#7 Lighted Small Circle Loop	626007
50 pk	#8 Lighted Narrow Spoon	626008
10 pk	#9 Lighted Tweezers	626009
50 pk	Variety Pack*	626010

*Six (6) each of 8 different curettes + 2 Tweezers

CURETTE MATERIAL POLYCARBONATE

INPUT 5V

INPUT CURRENT 150mA

WORKING CURRENT 10-15mA

POWER CONSUMPTION 0.05W

VOLTAGE 3.7-4.2V

EXECUTIVE STANDARD 4706.59-2008/GB 4706.1-2005

CERTIFICATION CE certificate based on the standard of EN55015

BATTERY CAPACITY 35mAh (approximately 50 minutes of use before needing to recharge)



POTENTIAL REIMBURSEMENT (U.S. ONLY)

CPT Code / Description

69210 – For potential reimbursement of impacted cerumen.



3796 Dunlap Street N, Arden Hills, MN 55112
1-800-328-4757 | info@warnertechcare.com
www.warnertechcare.com



LBEC SL5H 0122 Rev D