

Tel. (603) 899-3089 Fax. (603) 899-9747 www.grasonassociates.com

What eartips fit my instrument?



Please use this guide to find out which eartips fit your instrument.

The list is sorted by instrument make. Some instruments may be listed twice due to differences in probe or other optional features. Be aware that some instruments or makes may have a note associated with them.

To Search this Document:

Use the find function located under the edit menu or press ctrl-f. It is recommended to search by manufacturer and look at multiple results, not just the first. Some instruments fit multiple types of eartips depending on probe.



If unsure of which eartip fits, please ask for a sample to try.

Know of Other Instruments that single use eartips fit?

Please let us know at: sales@grasonassociates.com

Document Format will be as below:

Instrument Make

Instrument Name X · · · · · KR, RS, MF | Instrument Name Y · · · · · · IA

April 2017 Page 1 of 5

American Electromedics

Note: To our knowledge the KR, RS and MF series of eartips will fit the older American Electromedics Instruments, although fit is not exact. To Use KR, RS or MF eartips with an American Electromedics instrument, the eartips must be pushed onto the probe a long way until the forward edge of the probe is flush with the front of the eartip. The only known option for exact fit on these instruments is to find Maico Rock eartips.

| 81 · · · · · KR, RS, MF | 85 · · · · · KR, RS, MF |
|-------------------------|-------------------------|
| 83 ····· KR, RS, MF | 86 · · · · · KR, RS, MF |

Amplaid

Note: To our knowledge the IT series of eartips will fit Amplaid Instruments.

Amplivox

| Otowave models | Otowave models |
|---|-----------------------------|
| with White Probe · · · · · · KR, RS, MF | with Clear Probe · · · · IA |

Danplex

Note: To our knowledge the KR, RS and MF series of eartips will fit Danplex Instruments.

EAR

Note: To our knowledge the KR, RS and MF series of eartips will fit the EAR 3A Earphones with EAR3-06 adapter.

Echodia

| Elios · · · · · KR, RS, MF | Babyscreen ····· KR, RS, MF |
|-----------------------------|-----------------------------|
| Otodia · · · · · KR, RS, MF | |

Etymotic Research

| Ero Scan DPOAE ······ KR, RS, MF | ER10D ····· KR, RS, MF |
|---------------------------------------|------------------------------|
| ER10 Low Noise™ | ER3A Earphone with probe |
| Microphone ER10B · · · · · KR, RS, MF | tip adapter ····· KR, RS, MF |

Grason Stadler

| Tympstar····· KR, RS, MF | GSI 20 through 39····· KR, RS, MF |
|--------------------------|--|
| Tympstar Pro · · · · IA | 1720 through 1739 · · · · · KR, RS, MF |
| Audera ····· GS | Corti · · · · · IA |
| Audio Screener ···· GS | 70 · · · · · GS |

April 2017 Page 2 of 5

Intelligent Hearing Systems

Note: To our knowledge the KR, RS and MF series of eartips will fit the Smart OAE with 10D probe.

Interacoustics

Note: To our knowledge many Interacoustics instruments can fit two different eartips depending on the model vintage. As a general rule for AT235, AT235h, AA220, AA222, and MT10: clear probes use the IA eartip, while white, grey or black probes use the KR, RS, MF eartips. If unsure, please ask us or your distributor for a sample to check fit.

| Otoread | Eclipse OAE (NEOPT probe) IA Eclipse (clear probe) IA IMP440 IA MT10 (clear probe) IA MT10 KR, RS, MF MT5 IA MTP IA EP25 IA |
|--|---|
| <u>Inventis SRL</u> | |
| Flute ···· IT Clarinet ··· IT | Viola ····· IT |
| Note: To our knowledge the KR, RS and MF series of eartips will fit Kamplex Instruments. | |
| Labat Srl | |
| Echolab · · · · · KR, RS, MF | |
| Madsen | |
| Otoflex ····· MO Zodiac ···· IA Zodiac 901 ···· KR, RS, MF | Alpha PM Capella 2 KR, RS, MF Accuscreen IA |

April 2017 Page 3 of 5

Maico

Note: To our knowledge many Maico instruments can fit two different eartips depending on the model vintage. As a general rule for Otowave, EasyTymp, MI22, MI23, MI24, MI26, MI34, MI44: clear probes use the IA eartip, while white, grey or black probes use the KR, RS, MF eartips. If unsure, please ask us or your distributor for a sample to check fit.

| Otowave (clear probe) IA Otowave KR, RS, MF MA790 KR, RS, MF MA800 KR, RS, MF Eroscan KR, RS, MF Eroscan Pro KR, RS, MF Eroscan Pro + Tymp KR, RS, MF Eroscan Plus DPOAE IA MI 23 (clear probe) IA MI 23 KR, RS, MF MI 26 (clear probe) IA MI 26 KR, RS, MF 620 KR, RS, MF | EasyTymp (clear probe) IA EasyTymp KR, RS, MF 630 KR, RS, MF MI 22 (clear probe) IA MI 22 KR, RS, MF Race Car Tymp (clear probe) IA Race Car Tymp KR, RS, MF MI 24 (clear probe) IA MI 24 KR, RS, MF MI 34 (clear probe) IA MI 34 KR, RS, MF MI 34 (clear probe) IA MI 34 KR, RS, MF MI 44 (clear probe) IA MI 44 KR, RS, MF |
|--|--|
| MedRx Otomore (clear make) | |
| Otowave (clear probe) ······IA | |
| Micro Audiometrics | |
| Earscan Acoustic Impedance Audiometer | r······KR, RS, MF |
| Neurosoft | |
| Neuro-Audio-Screen ···· KR, RS, MF | Neuro-Audio/OAE · · · · · · KR, RS, MF |
| <u>Oscilla</u> | |
| TSM 300····· KR, RS, MF | TSM 400····· KR, RS, MF |

April 2017 Page 4 of 5

Otodynamics

Note: Grason & Associates sells two types of eartips that to our knowledge fit the Otodynamics instruments. TEOAE only probes that look like on the end will fit our OT0XX and OT2XX eartips. TEOAE + DPOAE probes that look like on the end will fit our OT3XX eartips.

| , | UGD probe ······ OT3XX SND probe ····· OT3XX |
|---|--|
| SNS probe ········· OT0XX, OT2XX SGS probe ····· OT0XX, OT2XX | SGD probe ············ OT3XX |

Path Medical

Resonance

Note: Eartip fit for the model numbers listed below will apply to all optional models. As an example, a listing for r25c will also fit the r25c-pt.

| r16m····· IT | r15c IT |
|--------------|-------------|
| r26m···· IT | r25c···· IT |
| r36m···· IT | r35c···· IT |

Siemens

Note: To our knowledge the KR, RS and MF series of eartips will fit Siemens Instruments.

<u>Virtual</u>

Note: To our knowledge the KR, RS and MF series of eartips will fit Virtual Instruments.

Vivosonic

Integrity OAE $\cdots \cdots$ KR, RS, MF

Welch Allyn

| 262 ····· KR, RS, MF | Micro Tymp ······WA |
|------------------------------------|-----------------------|
| 39500 OAE Screener ···· KR, RS, MF | Micro Tymp 2 ······WA |

April 2017 Page 5 of 5