

Otoscope Specula Guide

	Heine UniSpec	Heine AllSpec	Welch Allyn KleenSpec 524 Series	Welch Allyn Reusable 243 Series	Welch Allyn KleenSpec 521 Series	Welch Allyn Reusable 220 Series	Firefly Specula	Hotchkiss Specula
Heine K100	X							
Heine K180		X						
Heine Beta 400		X						
Heine Mini 3000		X						
LuxaScope		X						
Welch Allyn Diagnostic			X	X				
Welch Allyn PocketScope			X	X				
Welch Allyn MacroView			X	X				
Welch Allyn Pneumatic					X	X		
Welch Allyn Operating					X	X		
ADC Diagnostix			X	X				
Dino-Lite					X	X		
Firefly Video Otoscopes							X	
Hotchkiss								X
MedRx Video Otoscope		X						