Printing date 06.05.2019

Revision: 06.05.2019

SECTION 1: Identification of the substance/mixture and of the company/ undertaking

- · 1.1 Product identifier
- · Trade name: Otoform Kc
- **1.2 Relevant identified uses of the substance or mixture and uses advised against** No further relevant information available.
- · Application of the substance / the mixture Condensation-vulcanising ear impression silicone.
- · 1.3 Details of the supplier of the safety data sheet
- Manufacturer/Supplier: Dreve Otoplastik GmbH Max-Planck-Straße 31 59423 Unna / Germany Tel.: ++49 2303 / 8807-0 Fax.: ++49 2303 / 8807-29
- Further information obtainable from: Department Research & Development Fax: ++49 2303 / 8807-562 Email: sicherheitsdatenblatt@dreve.de
- **1.4 Emergency telephone number:** Tel.: ++ 49 211 / 797-3350 Plant Fire Department Henkel

SECTION 2: Hazards identification

- · 2.1 Classification of the substance or mixture
- Classification according to Regulation (EC) No 1272/2008 The product is not classified, according to the CLP regulation.

· 2.2 Label elements

- · Labelling according to Regulation (EC) No 1272/2008 Void
- · Hazard pictograms Void
- · Signal word Void
- · Hazard statements Void
- · Additional information:
- Safety data sheet available on request.
- · 2.3 Other hazards
- · Results of PBT and vPvB assessment
- · PBT: Not applicable.
- · vPvB: Not applicable.

SECTION 3: Composition/information on ingredients

- · 3.2 Chemical characterisation: Mixtures
- · Description: Condensation-vulcanising kneadable silicone

· Dangerous compo	nents:	
	White mineral oil, petroleum	10-25%
EINECS: 232-455-8	🕹 Asp. Tox. 1, H304	
	Paraffin waxes and Hydrocarbon waxes	≤2.5%
EINECS: 232-315-6	substance with a Community workplace exposure limit	
Additional information	tion: For the wording of the listed hazard phrases refer to section 16.	

(Contd. on page 2)

GR



Printing date 06.05.2019

Version number 3

Revision: 06.05.2019

Trade name: Otoform Kc

(Contd. of page 1)

SECTION 4: First aid measures

· 4.1 Description of first aid measures

- · General information: No special measures required.
- After inhalation: Supply fresh air; consult doctor in case of complaints.
- After skin contact: Generally the product does not irritate the skin.
- · After eye contact:

Rinse opened eye for several minutes under running water. Then consult a doctor.

- After swallowing: Do not induce vomiting. Call a doctor immediately.
- **4.2 Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed** No further relevant information available.

SECTION 5: Firefighting measures

- · 5.1 Extinguishing media
- Suitable extinguishing agents: CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- For safety reasons unsuitable extinguishing agents: Water with full jet
- 5.2 Special hazards arising from the substance or mixture No further relevant information available.
- 5.3 Advice for firefighters
- Protective equipment: No special measures required.

SECTION 6: Accidental release measures

- 6.1 Personal precautions, protective equipment and emergency procedures Avoid contact with eyes.
- · 6.2 Environmental precautions: Do not allow to enter sewers/ surface or ground water.
- 6.3 Methods and material for containment and cleaning up: Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- 6.4 Reference to other sections
 No dangerous substances are released.
 See Section 7 for information on safe handling.
 See Section 8 for information on personal protection equipment.
- See Section 13 for disposal information.

SECTION 7: Handling and storage

- **7.1 Precautions for safe handling** Avoid contact with eyes When using do not eat, drink or smoke.
- · Information about fire and explosion protection: No special measures required.
- · 7.2 Conditions for safe storage, including any incompatibilities
- · Storage:
- Requirements to be met by storerooms and receptacles: Keep container tightly sealed.
- · Information about storage in one common storage facility: Not required.
- Further information about storage conditions: None.

(Contd. on page 3)

⁻ GB



Printing date 06.05.2019

Version number 3

Revision: 06.05.2019

Trade name: Otoform Kc

· 7.3 Specific end use(s) No further relevant information available.

(Contd. of page 2)

SECTION 8: Exposure controls/personal protection

· Additional information about design of technical facilities: No further data; see item 7.

· 8.1 Control parameters

· Ingredients with limit values that require monitoring at the workplace:

8002-74-2 Paraffin waxes and Hydrocarbon waxes

WEL Short-term value: 6 mg/m³

Long-term value: 2 mg/m³

• Additional information: The lists valid during the making were used as basis.

· 8.2 Exposure controls

· Personal protective equipment:

- · General protective and hygienic measures:
- The usual precautionary measures are to be adhered to when handling chemicals.

• Respiratory protection: Not required.

• Protection of hands: Not required.

· Material of gloves Not applicable.

- · Penetration time of glove material Not applicable.
- Eye protection: Goggles recommended during refilling
- Body protection: Protective work clothing

SECTION 9: Physical and chemical properties · 9.1 Information on basic physical and chemical properties · General Information · Appearance: Form: Kneadable Colour: White · Odour: Characteristic · Odour threshold: Not determined. · pH-value: Not determined. · Change in condition Melting point/freezing point: Undetermined. Initial boiling point and boiling range: >300 °C >161 °C · Flash point: · Flammability (solid, gas): Not applicable. · Ignition temperature: >200 °C > 200 °C · Decomposition temperature: Auto-ignition temperature: Product is not selfigniting. · Explosive properties: Product does not present an explosion hazard. · Explosion limits: Lower: 0.6 Vol % 6.5 Vol % Upper: · Oxidising properties Not determined (Contd. on page 4)

GB



Printing date 06.05.2019

Version number 3

Revision: 06.05.2019

Trade name: Otoform Kc

	(Contd. of page
· Vapour pressure at 20 °C:	100 hPa
· Density at 20 °C:	1.4 g/cm ³
· Relative density	Not determined.
· Vapour density	Not determined.
· Evaporation rate	Not determined.
· Solubility in / Miscibility with	
water:	Not miscible or difficult to mix.
· Partition coefficient: n-octanol/water:	Not determined.
· Viscosity:	
Dynamic:	Not determined.
Kinematic:	Not determined.
· 9.2 Other information	No further relevant information available.

SECTION 10: Stability and reactivity

- · 10.1 Reactivity No dangerous reactions if used according to specifications.
- 10.2 Chemical stability Stable if used according to specifications.
- Thermal decomposition / conditions to be avoided:
- No decomposition if used according to specifications.
- 10.3 Possibility of hazardous reactions No dangerous reactions known.
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- · 10.6 Hazardous decomposition products: No dangerous decomposition products known.

SECTION 11: Toxicological information

- · 11.1 Information on toxicological effects
- Acute toxicity Based on available data, the classification criteria are not met.
- · Primary irritant effect:
- · Skin corrosion/irritation Based on available data, the classification criteria are not met.
- Serious eye damage/irritation Based on available data, the classification criteria are not met.
- · Respiratory or skin sensitisation Based on available data, the classification criteria are not met.
- · CMR effects (carcinogenity, mutagenicity and toxicity for reproduction)
- · Germ cell mutagenicity Based on available data, the classification criteria are not met.
- · Carcinogenicity Based on available data, the classification criteria are not met.
- · Reproductive toxicity Based on available data, the classification criteria are not met.
- · STOT-single exposure Based on available data, the classification criteria are not met.
- · STOT-repeated exposure Based on available data, the classification criteria are not met.
- Aspiration hazard Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

- · 12.1 Toxicity
- Aquatic toxicity: No further relevant information available.
- 12.2 Persistence and degradability No further relevant information available.
- · 12.3 Bioaccumulative potential No further relevant information available.
- 12.4 Mobility in soil No further relevant information available.

(Contd. on page 5)

GΒ



Printing date 06.05.2019

Version number 3

Revision: 06.05.2019

Trade name: Otoform Kc

(Contd. of page 4)

- · Additional ecological information:
- General notes: Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- · 12.5 Results of PBT and vPvB assessment
- **PBT:** Not applicable.
- · vPvB: Not applicable.
- **12.6 Other adverse effects** No further relevant information available.

SECTION 13: Disposal considerations

- · 13.1 Waste treatment methods
- Recommendation Disposal must be made according to official regulations.
- · Uncleaned packaging:
- Recommendation: Disposal must be made according to official regulations.

SECTION 14: Transport informati	on	
· 14.1 UN-Number · ADR, ADN, IMDG, IATA	Void	
 · 14.2 UN proper shipping name · ADR, ADN, IMDG, IATA 	Void	
· 14.3 Transport hazard class(es)		
· ADR, ADN, IMDG, IATA · Class	Void	
 14.4 Packing group ADR, IMDG, IATA 	Void	
 · 14.5 Environmental hazards: · Marine pollutant: 	No	
· 14.6 Special precautions for user	Not applicable.	
 14.7 Transport in bulk according to Anr of Marpol and the IBC Code 	nex II Not applicable.	
· UN "Model Regulation":	Void	

SECTION 15: Regulatory information

- \cdot 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture
- Directive 2012/18/EU
- Named dangerous substances ANNEX I None of the ingredients is listed.
- 15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

(Contd. on page 6)



Printing date 06.05.2019

Version number 3

Revision: 06.05.2019

Trade name: Otoform Kc

(Contd. of page 5)

he International Transport of Dangerous Goods by Rail) CAO: International Civil Aviation Organisation	earch & Development anne Langesberg
Department issuing SDS: Department Rese Contact: Dr. Thomas Veit, Lothar Sutor, Sus Abbreviations and acronyms: RID: Règlement international concernant le transport des m he International Transport of Dangerous Goods by Rail) CAO: International Civil Aviation Organisation	earch & Development
Contact: Dr. Thomas Veit, Lothar Sutor, Sus Abbreviations and acronyms: RID: Règlement international concernant le transport des m he International Transport of Dangerous Goods by Rail) CAO: International Civil Aviation Organisation	anne Langesberg
Contact: Dr. Thomas Veit, Lothar Sutor, Sus Abbreviations and acronyms: RID: Règlement international concernant le transport des m he International Transport of Dangerous Goods by Rail) CAO: International Civil Aviation Organisation	anne Langesberg
Abbreviations and acronyms: RID: Règlement international concernant le transport des m he International Transport of Dangerous Goods by Rail) CAO: International Civil Aviation Organisation	
RID: Règlement international concernant le transport des rr he International Transport of Dangerous Goods by Rail) CAO: International Civil Aviation Organisation	narchandises dangereuses par chemin de fer (Regulations Concern
he International Transport of Dangerous Goods by Rail) CAO: International Civil Aviation Organisation	
CAO: International Civil Aviation Organisation	
ADR: Accord europeen sur le transport des marchandis	ses dangereuses par Route (European Agreement concerning
nternational Carriage of Dangerous Goods by Road)	
MDG: International Maritime Code for Dangerous Goods	
ATA: International Air Transport Association	
GHS: Globally Harmonised System of Classification and La	
EINECS: European Inventory of Existing Commercial Chen	nical Substances
ELINCS: European List of Notified Chemical Substances	Chemical Society)
CAS: Chemical Abstracts Service (division of the American PBT: Persistent, Bioaccumulative and Toxic	
PvB: very Persistent and very Bioaccumulative	
Asp. Tox. 1: Aspiration hazard – Category 1	
⁷ Data compared to the previous version a	