

Printing date 12/10/2015

Reviewed on 12/10/2015

1 Identification

- · Product identifier
- · Trade name: Otoform pink siliclone
- · Application of the substance / the mixture Impression material
- · Details of the supplier of the safety data sheet
- · Manufacturer/Supplier:

Dreve Otoplastik GmbH Max-Planck-Straße 31 59423 Unna / Germany Tel.: ++49 2303 / 8807-0

Fax.: ++49 2303 / 8807-0

· Information department:

Department Research & Development

Fax.: +49 2303 8807-562

Email: sicherheitsdatenblatt@dreve.de • Emergency telephone number:

+49 211 797-3350

Plant Fire Department Henkel

2 Hazard(s) identification

· Classification of the substance or mixture

The product is not classified according to the Globally Harmonized System (GHS).

- · Label elements
- · GHS label elements Void
- · Hazard pictograms Void
- · Signal word Void
- · Hazard statements Void
- · Classification system:
- NFPA ratings (scale 0 4)



Health = 0 Fire = 1 Reactivity = 0

· HMIS-ratings (scale 0 - 4)



Health = 0 Fire = 1

Reactivity = 0

- · Other hazards
- · Results of PBT and vPvB assessment
- PBT: Not applicable.vPvB: Not applicable.

3 Composition/information on ingredients

- · Chemical characterization: Mixtures
- · Description: Addition-vulcanising 2-component silicone

· Dangerous components:		
14464-46-1	cristobalite	10-25%
471-34-1	calcium carbonate	2,5-10%

(Contd. on page 2)



Printing date 12/10/2015

Reviewed on 12/10/2015

Trade name: Otoform pink siliclone

	(Col	ntd. of page 1)
8042-47-5	paraffin oil	2,5-10%
14807-96-6	Talc (Mg3H2(SiO3)4)	≤ 2,5%
93685-79-1	Isoparaffin-mixture C12, C16 and C20	≤ 2,5%
8002-74-2	Paraffin waxes and Hydrocarbon waxes	≤ 2,5%

· Additional information: For the wording of the listed risk phrases refer to section 16.

4 First-aid measures

- · Description of first aid measures
- · General information: No special measures required.
- · After inhalation: Supply fresh air; consult doctor in case of complaints.
- · After skin contact: Generally the product does not irritate the skin.
- · After eye contact:

Rinse opened eye for several minutes under running water. Then consult a doctor.

· After swallowing:

Do not induce vomiting.

Call a doctor immediately.

- · Information for doctor:
- \cdot Most important symptoms and effects, both acute and delayed

No further relevant information available.

· Indication of any immediate medical attention and special treatment needed

No further relevant information available.

5 Fire-fighting measures

- · Extinguishing media
- · Suitable extinguishing agents:

CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

- · For safety reasons unsuitable extinguishing agents: Water with full jet
- · Special hazards arising from the substance or mixture

No further relevant information available.

- · Advice for firefighters
- · Protective equipment: No special measures required.

6 Accidental release measures

- Personal precautions, protective equipment and emergency procedures
 Avoid contact with eyes.
- Environmental precautions: Do not allow to enter sewers/ surface or ground water.
- · Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

· Reference to other sections

No dangerous substances are released.

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

- บร



Printing date 12/10/2015 Reviewed on 12/10/2015

Trade name: Otoform pink siliclone

(Contd. of page 2)

7 Handling and storage

- · Handling:
- **Precautions for safe handling** When using do not eat, drink or smoke.
- · Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
- · Storage:
- · Requirements to be met by storerooms and receptacles: Keep receptacle tightly sealed.
- · Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: None.
- · Specific end use(s) No further relevant information available.

8 Exposure controls/personal protection

- · Additional information about design of technical systems: No further data; see item 7.

· Cont	· Control parameters		
· Com	· Components with limit values that require monitoring at the workplace: 14464-46-1 cristobalite		
1446			
PEL	½ value from respirable dust formulae for Quartz		
REL	Long-term value: 0.05* mg/m³ *respirable dust; See Pocket Guide App. A		
TLV	Long-term value: 0.025* mg/m³ *as respirable fraction		
471-3	34-1 calcium carbonate		
PEL	Long-term value: 15* 5** mg/m³ *total dust **respirable fraction		
REL	Long-term value: 10* 5** mg/m³ *total dust **respirable fraction		
TLV	TLV withdrawn		
8002	8002-74-2 Paraffin waxes and Hydrocarbon waxes		
REL	Long-term value: 2 mg/m³		
TLV	Long-term value: 2 mg/m³		

- · Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment:
- · General protective and hygienic measures:

The usual precautionary measures for handling chemicals should be followed.

- · Breathing equipment: Not required.
- · Protection of hands: Not required.
- · Material of gloves Not applicable.
- · Penetration time of glove material Not applicable.
- · Eye protection: Goggles recommended during refilling.
- · Body protection: Protective work clothing



Printing date 12/10/2015 Reviewed on 12/10/2015

Trade name: Otoform pink siliclone

(Contd. of page 3)

nformation on basic physical and	chemical properties	
General Information		
Appearance:	\rangle ''	
Form:	Viscous	
Color: Odor:	Pink opaque Characteristic	
Odor threshold:	Not determined.	
pH-value:	Not determined.	
Change in condition		
Melting point/Melting range:	Undetermined.	
Boiling point/Boiling range:	301 °C (574 °F)	
Flash point:	131 °C (268 °F)	
Flammability (solid, gaseous):	Not applicable.	
Ignition temperature:	> 400 °C (> 752 °F)	
Decomposition temperature:	> 200 °C (> 392 °F)	
Auto igniting:	Product is not selfigniting.	
Danger of explosion:	Product does not present an explosion hazard.	
Explosion limits:		
Lower:	Not determined.	
Upper:	Not determined.	
Oxidizing properties	Not determined	
Vapor pressure at 20 °C (68 °F):	5 hPa (4 mm Hg)	
Density at 20 °C (68 °F):	1,5 g/cm³ (12.518 lbs/gal)	
Relative density	Not determined.	
Vapor density	Not determined.	
Evaporation rate	Not determined.	
Solubility in / Miscibility with		
Water:	Not miscible or difficult to mix.	
· Partition coefficient (n-octanol/water): Not determined.		
Viscosity:		
Dynamic:	Not determined.	
Kinematic:	Not determined.	
Solvent content:		
Organic solvents:	0,0 %	
Other information	No further relevant information available.	

10 Stability and reactivity

- · Reactivity No dangerous reactions if used according to specifications.
- · Chemical stability Stable if used according to specifications.
- Thermal decomposition / conditions to be avoided:
 No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known.

(Contd. on page 5)



Printing date 12/10/2015 Reviewed on 12/10/2015

Trade name: Otoform pink siliclone

(Contd. of page 4)

- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- · Information on toxicological effects
- · Acute toxicity:
- · Primary irritant effect:
- · on the skin: No irritant effect.
- · on the eye: No irritant effect.
- · Sensitization: No sensitizing effects known.
- · Additional toxicological information:

The product is not subject to classification according to internally approved calculation methods for preparations:

When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

· Carcinogenic categories

· IARC (International Agency for Research on Cancer)		
14464-46-1	cristobalite	1
14807-96-6	Talc (Mg3H2(SiO3)4)	3
NTP (National Toxicology Program)		
14464-46-1	cristobalite	K
· OSHA-Ca (Occupational Safety & Health Administration)		
None of the ingredients is listed.		

12 Ecological information

- Toxicity
- · Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behavior in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes:

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- · Results of PBT and vPvB assessment
- · **PBT**: Not applicable.
- · **vPvB**: Not applicable.
- · Other adverse effects No further relevant information available.

13 Disposal considerations

- · Waste treatment methods
- · Recommendation: Disposal must be made according to official regulations.

(Contd. on page 6)



Printing date 12/10/2015 Reviewed on 12/10/2015

Trade name: Otoform pink siliclone

(Contd. of page 5)

- · Uncleaned packagings:
- · Recommendation: Disposal must be made according to official regulations.

14 Transport information	
· UN-Number · DOT, ADR, ADN, IMDG, IATA	Void
UN proper shipping nameDOT, ADR, ADN, IMDG, IATA	Void
· Transport hazard class(es)	
· DOT, ADR, ADN, IMDG, IATA · Class	Void
Packing group DOT, ADR, IMDG, IATA	Void
Environmental hazards: Marine pollutant:	No
· Special precautions for user	Not applicable.
· Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code Not applicable.	
· UN "Model Regulation":	Void

15 Regulatory information

- · Safety, health and environmental regulations/legislation specific for the substance or mixture
- · Sara

· Sara	
· Section 355 (extremely hazardous su	bstances):

None of the ingredient is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

· TSCA (Toxic Substances Control Act):

All ingredients are listed.

- · Proposition 65
- · Chemicals known to cause cancer:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

(Contd. on page 7)



Printing date 12/10/2015 Reviewed on 12/10/2015

Trade name: Otoform pink siliclone

(Contd. of page 6)

· Cancerogenity categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

· TLV (Threshold Limit Value established by ACGIH)

14807-96-6 Talc (Mg3H2(SiO3)4)

A4

· MAK (German Maximum Workplace Concentration)

14464-46-1 cristobalite

1

· NIOSH-Ca (National Institute for Occupational Safety and Health)

14464-46-1 cristobalite

- · GHS label elements Void
- · Hazard pictograms Void
- · Signal word Void
- Hazard statements Void
- · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- · Department issuing SDS: Department Research & Development
- · Contact: Dr. Thomas Veit, Lothar Sutor
- · Date of preparation / last revision 12/10/2015 / 1
- · Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

US