## Otoform® Ak X



## **Technical Data Sheet**

Intended use: Addition vulcanizing ear impression silicone

Manufacturer: Dreve Otoplastik GmbH

Max-Planck-Straße 31 59423 Unna / Germany

Composition:

Pos	chemical term
1	vinyl containing polydimethylsiloxane
2	silicic acid
3	filler materials
4	pigments
5	platinum catalyst

Color: white / pink

Processing time:  $100 \pm 15 \text{ sec}$ 

Setting time: 3 min  $\pm$  15 sec

Final hardness:  $30 \pm 2$  Shore A

Consistency: putty consistency - type 0

(according to DIN EN ISO 4823)

Recovery from deformation: >99,6 % (according to DIN EN ISO 4823)

Strain in compression: 4,4 % (according to DIN EN ISO 4823)

Linear dimensional change: <0,2 % (according to DIN EN ISO 4823)

Dreve Otoplastik GmbH

Max-Planck-Str. 31

59423 Unna / Germany

Tel.: 02303 / 8807-0

Fax: 02303 / 82909

otoplastik@dreve.de

www.dreve.com