

Hardener Paste



Technical Data Sheet

Intended use: Hardener paste for condensation vulcanizing impression materials for ear impression taking

Manufacturer: Dreve Otoplastik GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	chemical term
1	organic tin compounds
2	silicic acid ester
3	vinyltrimethoxysilane
4	food colour
5	silicon dioxide
6	vaseline

Form: paste

Colour: red
(corresponds to released sample in the house)

Density: 1.02 g/cm³ at 20 °C DIN 51757

Solubility in water: insoluble

pH-value: neutral

Dreve Otoplastik GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 82909
otoplastik@dreve.de
www.dreve.com