

# MICROSONIC, INC.

2960 DUSS AVENUE, AMBRIDGE, PENNSYLVANIA 15003 Tel: (800) 523-9270 Fax: (724) 266-6309

## MATERIAL SAFETY DATA SHEET BLUE - Mega Sil/Silicone Impression Material

### Section I - Product Information:

**Chemical Name:** Silicone Impression Material

**DOT/UN Shipping Name:** Vinyl Polysiloxane

**Publication Date:** January, 2012

**CAS Reg. No.:**

**Trade Name/Product Code:**

**Update:** January, 2019

**Formulation:** Proprietary Formulation with vinyl polysiloxanes. Catalyst with platinum complex and base with hydrogen vinyl polysiloxane. Inert fillers in both catalyst and base.

### Section II - Hazardous Ingredients:

Hazardous Components/Chemical Identity	CAS Number	%ACGIHTLV	OSHA PEL
NONE	NA		

### Section III-Physical Characteristics:

Boiling Point	Specific Gravity	Vapor Pressure	Percent Volatile	Vapor Density	Evaporation Rate	Solubility in Water
Non Volatile	(H <sub>2</sub> O=1): 1.55	NA	w/w%: Negligible	( Air=1): NA	(H <sub>2</sub> O): Negligible	Negligible

**Appearance and Odor:** Base has Blue Shade and no odor. The catalyst is Colorless and odorless.

### Section IV - Fire and Explosion Data:

Flash Point	Flammable Limit	Auto-ignition Temperature
NA	Air vol. %: Lower: NA; Upper: NA	NA
<b>Extinguisher Method:</b>	Water, Carbon Dioxide, Dry Chemicals	

**Special Fire Fighting Procedures:** None

**Explosion Hazard:** None

**Unusual Fire and Explosion Hazards:** None

### Section V: Composition / Information on Ingredients

<u>Chemical Identity</u>	<u>% Concentration</u>
Silicone Polymers - Fully and/or partially crosslinked	>83%
Fillers - Inert	>10%
Minerals & Pigments	<2%
Solvent - Mineral oils	>5%

### Section VI- Reactivity:

**Stability:** Stable

**Conditions to Avoid:** Do not use with sulfur containing latex gloves. Do not use with materials containing amines

**Incompatibility (Materials to Avoid):** Strong oxidizing agents

**Hazardous Decomposition Products:** None

**Hazardous Polymerization:** Will not occur

**Conditions to Avoid:** None

**Section VII-Health Hazard data:**

**Primary Routes of Entry:** Skin and Eyes

**Carcinogenicity:** None of the components of this material are listed by IARC, NTP, OSHA, or ACIGH as carcinogens.

**Target Organs:** For the mixture: NA.

**Threshold Limit Value(LTV):** For the mixture: NA

**Effects of Over Exposure:** It is not known to cause significant health problems.

**Toxicity Information:** For the mixture: None Listed. All ingredients are listed on the TOSCA inventory.

**Section VI - Health Hazard data: (Continued):**

**Emergency and First Aid Procedure:**

- Inhalation:** Remove to fresh air. If breathing has stopped, give oxygen. Get medical help if discomfort persists.
- Eyes:** Flush with water for 15 minutes, including under eyelids. Get medical help if discomfort persists.
- Skin:** Wash thoroughly with soap and water. Remove contaminated clothing. Get medical help if discomfort persists.

**Ingestion:** Rinse mouth out with water. Call a physician if the amount was large.

**Hazardous Material Identification System (HMIS) Rating:**

**Health:** 0 **Flammability:** 0 **Reactivity:** 0 **Personal Protective Equipment:** Gloves, splash goggles

**Section VII - Spill or Leak Procedure:**

**Step to Be Taken in Case Material is Released or Spilled:** Dike and sweep up to avoid slipping hazards.

**Waste Disposal Method:** may be disposed of in a landfill or incinerated. Follow Federal, State, and Local regulations for disposal.

**Section VIII - Precautions for Safe Handling and Use:**

**Respiratory Protection:** None required.

**Ventilation:** A good exhaust fan.

**Protective Gloves:** Latex gloves which contain sulfur will inhibit polymerization of two component systems. Avoid the use of this style of a glove.

**Eye Protection:** Safety Glasses or Splash Goggles.

**Other Protective Equipment:** Provide an eyewash and safety shower.

**Section IX - Special Precaution:**

**Precaution to Be Taken in Handling and Storage:** Observe precautions found on the label. Store in cool dry place away from heat, sparks, flame and direct sunlight. Close the container after each use. Ground all metal containers when transferring.

**Other Precautions:** Wash face and hands thoroughly with soap and water after handling and before eating, drinking or smoking

=====

The information presented herein was obtained from sources considered to be reliable. However, this information is provided without any warranty, expressed or implied, regarding its correctness or suitability for consumers intended use and/or application. For this and other reasons, we assume no responsibility and expressly disclaim liability for loss, damage or expense arising out of or in any way connected with the handling, storage, use or disposal of the product. This MSDS was prepared expressly for this product. Use the materials only as directed. If the product is used as a component of another product, the information contained within the MSDS may not be applicable.

