

Product Name: Cavicide

Product Description Form: Hard surface cleaner and disinfectant

Application: All fill sizes, see www.metrex.com for catalogue numbers and safety data sheets.

Company: METREX® RESEARCH
1717 W. Collins Ave.
Orange, CA 92867
U.S.A.
1-800-841-1428

Date of Disclosure

February 16, 2020

Ingredients

Ingredient Name	CAS	Present on Designated List	Function (HCPA)
Water	7732-18-5	No	Solvent
Isopropanol	67-63-0	No	Solvent / active ingredient
Ethylene Glycol Monobutyl Ether	111-76-2	CA TACS	Solvent
Diisobutylphenoxyethoxyethyl dimethyl benzyl ammonium chloride	121-54-0	No	Active ingredient
Sodium Lauryl Ether Sulfate	9004-82-4	No	Cleaner
Surfactant	With held for CBI	EU Endocrine Disruptors, Marine Priority Action Chemical	Cleaner
Sodium Hydroxide	1310-73-2	No	pH adjustment

Contact for Additional Ingredient Information

Links

CA TACS <https://ww2.arb.ca.gov/resources/documents/carb-identified-toxic-air-contaminants>

EU Endocrine Disruptors https://echa.europa.eu/candidate-list-table?p_p_id=disslists_WAR_disslistsportlet&p_p_lifecycle=0&p_p_state=normal&p_p_mode=view&p_p_col_id=column-1&p_p_col_pos=2&p_p_col_count=3&disslists_WAR_disslistsportlet_keywords=&disslists_WAR_disslistsportlet_o

Marine Priority Action Chemical <https://www.ospar.org/work-areas/hasec/chemicals/priority-action>