



Product Name: egger LP/H modelling paste, red, blue, transparent

REF: 29610, 29615, 30800

Date of print: 01.02.2022 Date of last alteration: 01.02.2022

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

egger LP/H modelling paste - red, blue, transparent

| | |
|-------------|--------------------------|
| red | UFI: GT8V-Y91Y-Q007-CFQU |
| blue | UFI: 2W8V-F9RD-000Q-1TA1 |
| transparent | UFI: QKDV-79EF-6000-GGCV |

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture

Ligth curing modeling paste for use in audiology.

1.3. Details of the supplier of the safety data sheet

Manufacturer/distributor: egger Otoplastik+Labortechnik GmbH

Street/POB-No.: Aybühlweg 59

Postal code/city/country: 87439 Kempten/Germany

Telephone: +49 831 58113-60

Telefax: +49 831 58113-14

Internet: www.egger.online

E-mail: sales@egger.online

1.4. Emergency telephone number: +49 89 19240

Emergency number: Munich, Germany, (toxicological dep. of the II. Med. Hospital)

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Regulation (EC) No. 1272/2008

Hazard categories:

Respiratory or skin sensitisation: Skin Sens. 1

Hazardous to the aquatic environment: Aquatic Chronic 2

Hazard Statements:

May cause an allergic skin reaction.

Toxic to aquatic life with long lasting effects.

2.2. Label elements

Regulation (EC) No. 1272/2008

Hazard components for labelling

7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxa-5,12-diazahexadecane-1,16-diyl bismethacrylate

diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide

Signal word: Warning

Pictograms:



Hazard statements

H317 May cause an allergic skin reaction.

H411 Toxic to aquatic life with long lasting effects.



Product Name: egger LP/H modelling paste, red, blue, transparent

REF: 29610, 29615, 30800

Date of print: 01.02.2022 Date of last alteration: 01.02.2022

Precautionary statements

- P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
 P273 Avoid release to the environment.
 P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.
 P302+P352 IF ON SKIN: Wash with plenty of water.
 P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
 P362+P364 Take off contaminated clothing and wash it before reuse.
 P391 Collect spillage.
 P501 Dispose of contents/container to according to local and applicable legislation of dispose of waste.

Special labelling of certain mixtures

EUH210 Safety data sheet available on request.

2.3. Other hazards

No information available.

SECTION 3: Composition/information on ingredients

3.2. Mixtures

Chemical characterization

Mixture of acrylic/methacrylic resins with auxilliary matters.

Hazardous components

| CAS No | Chemical name | | | Quantity |
|------------|--|--------------|------------------|-------------|
| | EC no | Index No | REACH No | |
| | GHS Classification | | | |
| 72869-86-4 | 7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxo-5,12-diazahexadecane-1,16-diyl bismethacrylate | | | 30 - < 75 % |
| | 276-957-5 | | 01-2120751202-68 | |
| | Skin Sens. 1B, Aquatic Chronic 2; H317 H411 | | | |
| 75980-60-8 | diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide | | | < 1 % |
| | 278-355-8 | 015-203-00-X | | |
| | Repr. 2, Skin Sens. 1B, Aquatic Chronic 2; H361f H317 H411 | | | |

Full text of H and EUH statements: see section 16.

SECTION 4: First aid measures

4.1. Description of first aid measures

After inhalation

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

After contact with skin

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water. In case of eye irritation consult an ophthalmologist.

After ingestion

Rinse mouth immediately and drink plenty of water.



**Product Name: egger LP/H modelling paste,
red, blue, transparent**

REF: 29610, 29615, 30800

Date of print: 01.02.2022 Date of last alteration: 01.02.2022

egger



Seek immediately medical advice. Do not induce vomiting. In case of spontaneous vomiting take care of an unhindered flow out of the vomit (danger of suffocation).

4.2. Most important symptoms and effects, both acute and delayed

No information available.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

5.2. Special hazards arising from the substance or mixture

Non-flammable.

5.3. Advice for firefighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately.

Do not allow entering drains or surface water.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

No special measures are necessary.



Product Name: egger LP/H modelling paste, red, blue, transparent

REF: 29610, 29615, 30800

Date of print: 01.02.2022 Date of last alteration: 01.02.2022

egger

Advice on protection against fire and explosion

No special fire protection measures are necessary.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep container tightly closed. Keep only in the original container in a cool, dry and well-ventilated place, away from foodstuffs. Keep away from all kind of lighth. An inert gas blanket should not be applied, because the stability of the product depends on the presence of oxygen (air).

Hints on joint storage

Keep away from spontaneous flammable or combustible substances.

Further information on storage conditions

Keep away from all kind of lighth. An inert gas blanket should not be applied, because the stability of the product depends on the presence of oxygen (air).

7.3. Specific end use(s)

Lighth curing material for repair and modeling of earmoulds.

For use by trained specialist staff.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

8.2. Exposure controls

Protective and hygiene measures

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary . When using do not eat, drink, smoke, sniff.

Eye/face protection

Wear eye/face protection.

Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

Suitable are gloves of the following material: Butyl caoutchouc (butyl rubber).

Skin protection

Wear suitable protective clothing.

Respiratory protection

In case of inadequate ventilation wear respiratory protection.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state: paste , low-viscosity

Colour: several shades, according to the product name

Odour: faintly like esters



Product Name: egger LP/H modelling paste, red, blue, transparent

REF: 29610, 29615, 30800

Date of print: 01.02.2022 Date of last alteration: 01.02.2022



| | | Test method |
|--|--------------------------------|-------------|
| pH-Value: | not determined | |
| Changes in the physical state | | |
| Melting point: | not determined | |
| Initial boiling point and boiling range: | not determined | |
| Flash point: | >100 °C | DIN 51755 |
| Flammability | | |
| Solid: | not applicable | |
| Gas: | not applicable | |
| Explosive properties | The product is not: Explosive. | |
| Lower explosion limits: | not determined | |
| Upper explosion limits: | not determined | |
| Auto-ignition temperature | | |
| Solid: | not applicable | |
| Gas: | not applicable | |
| Decomposition temperature: | >=190 °C | |
| Oxidizing properties | Not oxidizing. | |
| Vapour pressure (at 20 °C): | <1 hPa | |
| Density (at 20 °C): | 1.25-1.35 g/cm ³ | DIN 51757 |
| Water solubility: | insoluble | |
| Solubility in other solvents | not determined | |
| Partition coefficient: | not determined | |
| Vapour density: | not determined | |
| Evaporation rate: | not determined | |

9.2. Other information

Solid content: not determined

SECTION 10: Stability and reactivity

10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

10.3. Possibility of hazardous reactions

Reacts with: strong oxidising agents, strong alkaline or acidic materials.

10.4. Conditions to avoid

Ultra-violet light and daylight initiate polymerisation of the product. Therefore keep only in tightly closed containers away from any sources of light. Keep in a refrigerator at 2°C - 25°C / 36°F - 77 °F.

10.5. Incompatible materials

No information available.

10.6. Hazardous decomposition products

In case of fire, acrid acrylic fumes may occur.



Product Name: egger LP/H modelling paste, red, blue, transparent

REF: 29610, 29615, 30800

Date of print: 01.02.2022 Date of last alteration: 01.02.2022

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity

Based on available data, the classification criteria are not met.

For the product itself no toxicological data are available. In products with a comparable composition, a LD50 (orally, species rat) of > 5000 mg/kg has been found.

| CAS No | Chemical name | | | | |
|-------------------|---|------------------|---------|----------|--------|
| | Exposure route | Dose | Species | Source | Method |
| 72869-86-4 | 7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxa-5,12-diazahexadecane-1,16-diyl bismethacrylate | | | | |
| | oral | LD50 >5000 mg/kg | Rat | OECD 401 | |
| | dermal | LD50 >2000 mg/kg | Rat | OECD 402 | |
| 75980-60-8 | diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide | | | | |
| | oral | LD50 >5000 mg/kg | Rat | | |
| | dermal | LD50 >2000 mg/kg | Rat | | |

Irritation and corrosivity

Based on available data, the classification criteria are not met.

Sensitising effects

May cause an allergic skin reaction.

(7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxa-5,12-diazahexadecane-1,16-diyl bismethacrylate;
diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide)

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

STOT-single exposure

Based on available data, the classification criteria are not met.

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Additional information on tests

This mixture is classified as not hazardous according to regulation (EC) 1272/2008 [CLP].

SECTION 12: Ecological information

12.1. Toxicity

Toxic to aquatic life with long lasting effects.

| CAS No | Chemical name | | | | | |
|-------------------|---|----------------|-----------|-----------------------------------|----------|--------|
| | Aquatic toxicity | Dose | [h] [d] | Species | Source | Method |
| 72869-86-4 | 7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxa-5,12-diazahexadecane-1,16-diyl bismethacrylate | | | | | |
| | Acute crustacea toxicity | EC50 >1.2 mg/l | 48 h | Daphnia magna (Big water flea) | OECD 202 | |



Product Name: egger LP/H modelling paste, red, blue, transparent

REF: 29610, 29615, 30800

Date of print: 01.02.2022 Date of last alteration: 01.02.2022

| 75980-60-8 | diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide | | | | |
|------------|---|------------------|------|--------------------------------|--|
| | Acute algae toxicity | ErC50 >2.01 mg/l | 72 h | Scenedesmus subspicatus | |
| | Acute crustacea toxicity | EC50 3.53 mg/l | 48 h | Daphnia magna (Big water flea) | |
| | Acute bacteria toxicity | (>1000 mg/l) | 3 h | Activated sludge | |

12.2. Persistence and degradability

The product has not been tested.

| CAS No | Chemical name | | | |
|------------|---|-------|----|--------|
| | Method | Value | d | Source |
| | Evaluation | | | |
| 75980-60-8 | diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide | | | |
| | | 0-10% | 28 | |
| | Not readily biodegradable (according to OECD criteria). | | | |

12.3. Bioaccumulative potential

The product has not been tested.

Partition coefficient n-octanol/water

| CAS No | Chemical name | Log Pow |
|------------|---|---------|
| 75980-60-8 | diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide | 3.1 |

BCF

| CAS No | Chemical name | BCF | Species | Source |
|------------|---|-------|-------------------------------|--------|
| 75980-60-8 | diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide | 47-55 | Cyprinus carpio (Common Carp) | |

12.4. Mobility in soil

The product has not been tested.

12.5. Results of PBT and vPvB assessment

Not identified as PBT/vPvB substances.

12.6. Other adverse effects

No information available.

Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal recommendations

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

Contaminated packaging

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

**Product Name: egger LP/H modelling paste,
red, blue, transparent**

REF: 29610, 29615, 30800

Date of print: 01.02.2022 Date of last alteration: 01.02.2022

SECTION 14: Transport information**Land transport (ADR/RID)**

14.1. UN number: UN 3082
14.2. UN proper shipping name: ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S.
14.3. Transport hazard class(es): 9
14.4. Packing group: III
Hazard label: 9
Classification code: M6
Special Provisions: 274 335 375 601
Limited quantity: 5 L
Excepted quantity: E1
Transport category: 3
Hazard No: 90
Tunnel restriction code: -

Inland waterways transport (ADN)

14.1. UN number: UN 3082
14.2. UN proper shipping name: ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S.
14.3. Transport hazard class(es): 9
14.4. Packing group: III
Hazard label: 9
Classification code: M6
Special Provisions: 274 335 375 601
Limited quantity: 5 L
Excepted quantity: E1

Marine transport (IMDG)

14.1. UN number: UN 3082
14.2. UN proper shipping name: ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S.
14.3. Transport hazard class(es): 9
14.4. Packing group: III
Hazard label: 9
Classification code: M6
Special Provisions: 274, 335, 969
Limited quantity: 5 L
Excepted quantity: E1
EmS: F-A, S-F

Air transport (ICAO-TI/IATA-DGR)

14.1. UN number: UN 3082
14.2. UN proper shipping name: ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S.
14.3. Transport hazard class(es): 9
14.4. Packing group: III
Hazard label: 9
Special Provisions: A97 A158 A197
Limited quantity Passenger: 30 kg G
Passenger LQ: Y964
Excepted quantity: E1
IATA-packing instructions - Passenger: 964
IATA-max. quantity - Passenger: 450 L
IATA-packing instructions - Cargo: 964
IATA-max. quantity - Cargo: 450 L



Product Name: egger LP/H modelling paste, red, blue, transparent

REF: 29610, 29615, 30800

Date of print: 01.02.2022 Date of last alteration: 01.02.2022

14.6. Special precautions for user

No dangerous good in sense of this transport regulation.

14.7. Transport in bulk according to Annex II of MARPOL and the IBC Code

No dangerous good in sense of this transport regulation.

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU regulatory information

Information according to 2012/18/EU
(SEVESO III):

E2 Hazardous to the Aquatic Environment

National regulatory information

Employment restrictions: Observe restrictions to employment for juveniles according to the 'juvenile work protection guideline' (94/33/EC).

Water hazard class (D): 1 - slightly hazardous to water

Skin resorption/Sensitization: Causes allergic hypersensitivity reactions.

15.2. Chemical safety assessment

Chemical safety assessments for substances in this mixture were not carried out.

SECTION 16: Other information

Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route
(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

CLP: Classification, labelling and Packaging

REACH: Registration, Evaluation and Authorization of Chemicals

GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals

UN: United Nations

DNEL: Derived No Effect Level

DMEL: Derived Minimal Effect Level

PNEC: Predicted No Effect Concentration

ATE: Acute toxicity estimate

LL50: Lethal loading, 50 %

EL50: Effect loading, 50 %

EC50: Effective Concentration 50 %

ErC50: Effective Concentration 50 %, growth rate

NOEC: No Observed Effect Concentration

BCF: Bio-concentration factor

PBT: persistent, bioaccumulative, toxic

vPvB: very persistent, very bioaccumulative



Product Name: egger LP/H modelling paste, red, blue, transparent

REF: 29610, 29615, 30800

Date of print: 01.02.2022 Date of last alteration: 01.02.2022



RID: Regulations concerning the international carriage of dangerous goods by rail
 ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
 (Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)
 EmS: Emergency Schedules
 MFAG: Medical First Aid Guide
 ICAO: International Civil Aviation Organization
 MARPOL: International Convention for the Prevention of Pollution from Ships
 IBC: Intermediate Bulk Container
 SVHC: Substance of Very High Concern
 @1602.B016012

Classification for mixtures and used evaluation method according to Regulation (EC) No. 1272/2008 [CLP]

| Classification | Classification procedure |
|-------------------------|--------------------------|
| Skin Sens. 1; H317 | Calculation method |
| Aquatic Chronic 2; H411 | Calculation method |

Relevant H and EUH statements (number and full text)

H317 May cause an allergic skin reaction.
 H361f Suspected of damaging fertility.
 H411 Toxic to aquatic life with long lasting effects.
 EUH210 Safety data sheet available on request.

Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)