

## PRODUCT COMPARISON









dry-cap uv®

for end consumers

dry-cap uv 2® for end consumers

dry-cap uv 3® for end consumers dry-cap pro® for professionals

dry-cap uv® for end consumers

Operating mode:

Dimensions: (L x W x H) in mm

Short-time operation 1h Short-time operation 2h Short-time operation 2h for domestic use 105 x 105 x 93

for domestic use

for domestic use 144 x 82 x 96

Continuous operation for professionals

Short-time operation 1h for domestic use

Fast drying by efficient

120 x 120 x 120

90 x 90 x 38



air flow blower

Quality fan, ball bearing, 30.000 h, 45 dB(A), 7.000 rpm, Airflow 2,56 CFM

ality fan, ball bearing, 30.000 h, 36 dB(A), 7.000 rpm, Airflow 2,56 CFM

uality fan, ball bearing, 30.000 h, 18.5 dB(A), 7.000 rpm, Airflow 2,88 CFM

High-class fan, ball bearing, 50.000 h, 42 dB(A), 7.000 rpm, Airflow 2,74 CFM

105 x 105 x 93 + Hygiene-Pad

Quality fan, ball bearing, 30.000 h, 43 dB(A), 7.000 rpm, Airflow 2,56 CFM



Hygienic cleansing by UV-C light for your health

ıration: 1 min at start + 1 min at end lifetime: 15.000 h

ration: 2 min at start + 2 min at end lifetime: 15.000 h

ration: 2 min at start + 2 min at end lifetime: 15.000 h

ration: 1 min at start + 1 min at end lifetime: 30.000 h





Health improving blue light













Easy handling by simple touch on tactile buttons













Mobility by 2 x 5V USB standard



**Dual-USB** power supply

LiPo battery

**Dual-USB** power supply

USB power supply



Gravity sensor for safety shut-off













Automatic switch-off at cycle end

Switch-off after 30 min Switch-off after 2 h

Switch-off after 2 h

Switch-off after 30 min and acoustic signal

Switch-off after 30 min



All features controlled by microprocessor













Highest product quality with 2 year warranty













Energy-efficient gentle and absolutely safe









